PEDIATRIC CHEST 3T

INSTRUCTIONS AND PATIENT PREPARATION

Instructions	No NPO requirement. Free-breathing. Cover lung apices through mid-neck.
Position	Supine, head first
Coil	dStream Torso coil
Landmark	Center at mid chest

CONTRAST AGENT PARAMETERS

Contrast Agent	Dose	Injection Rate	Dose Timing Method and Scan Delay
Gadopiclenol (VUEWAY)	0.1 mL/kg (single dose, weight-based)	Hand injection	IV placed prior to start of MRI
Saline Flush	5-10 mL	Hand injection	Following contrast
Comments	Single dose of VUEWAY 0.1 mL/kg based on patient weight; hand inject followed by 5-10 mL saline flush. Report agent administered and dose.		

Disclaimer: The MR Protocols information is provided by independent physicians and technologists and is not intended to provide definitive or complete medical instructional or treatment information. Products referenced herein should be used consistently with their package inserts, and any devices referenced herein should be used consistently with the operation instructions provided by their respective manufactures. Northwest Imaging Forums, Inc. and Bracco Diagnostics Inc. do not offer medical advice and shall not be responsible for any person's reliance on this protocol.

Supported by an educational grant from Bracco Diagnostics Inc.

MR Gadopiclenol (VUEWAY®) Protocol for Philips MRI Scanners
Pediatric

NW IMAGING FORUMS, INC.

INFO PAGE

Scan Description	Survey	AX T2 FS MVXD	CO T2 FS MVXD	SG T2 FS MVXD	AX T1 FS mDIXON VANE 3D	AX T1 FS mDIXON VANE 3D POST CONTRAST	CO T1 FS 3D POST CONTRAST
Total scan duration	00:21.3	03:15.0	02:06.0	02:48.0	02:30.5	02:30.5	02:01.3
Rel. SNR	6.40643835	0.289420217	0.603155911	0.791086435	0.76686883	0.76686883	0.67733866
Act. TR(ms)	—	2725	1860	1860	—	—	—
Act. TE(ms)	—	149 (146/151)	80	80	—	—	—
Act. TR/TE(ms)	8.0 / 4.6	—	—	—	—	—	—
Act. TR/TI(ms)	—	—	—	—	—	—	—
Act. TR/TE1/delta TE(ms)	—	—	—	—	—	—	—
ACQ matrix M x P	172 x 85	184 x 184	156 x 156	156 x 156	140 x 140	140 x 140	140 x 140
ACQ voxel MPS (mm)	1.74 / 3.53 / 15.0	1.09 / 1.09 / 4.00	1.28 / 1.28 / 3.00	1.28 / 1.28 / 4.00	1.29 / 1.29 / 1.30	1.29 / 1.29 / 1.30	1.29 / 1.29 / 1.30
REC voxel MPS (mm)	1.56 / 1.56 / 15.0	0.83 / 0.83 / 4.00	0.89 / 0.89 / 3.00	0.89 / 0.89 / 4.00	1.02 / 1.02 / 1.00	1.02 / 1.02 / 1.00	1.02 / 1.02 / 1.00
Scan percentage (%)	49.4186058	—	—	—	—	—	—
Packages	—	5	6	7	—	—	—
TFE shots	—	—	—	—	—	—	—
TFE dur. shot / acq (ms)	708.8 / 677.0	—	—	—	663.1 / 564.2	663.1 / 564.2	532.2 / 438.6
TFE shot interval (ms)	708.826477	—	—	—	—	—	—
Min. TI delay	—	—	—	—	—	—	—
Min. slice gap (mm)	—	4	3	4	—	—	—
EPIfactor	—	—	—	—	—	—	—
Act.slicegap (mm)	—	—	—	—	-0.299999952	-0.299999952	-0.301075339
Act. WFS (pix) / BW (Hz)	2.107 / 206.2	—	—	—	0.452 / 960.1	0.452 / 960.1	0.452 / 960.1
Min. WFS (pix) / Max. BW (Hz)	0.377 / 1153.5	—	—	—	0.380 / 1144.7	0.380 / 1144.7	0.394 / 1102.3
Min. TR/TE(ms)	8.0 / 3.5	—	—	—	—	—	—
Min. TR/TE1/delta TE(ms)	—	—	—	—	—	—	—
WFS (pix) / BW (Hz)	—	0.422 / 1029.2	0.358 / 1214.0	0.358 / 1214.0	—	—	—
TSE es / shot(ms)	—	5.8 / 297	5.5 / 155	5.5 / 155	—	—	—
TEeff / TEequiv(ms)	—	146 / 97	80 / 58	80 / 58	—	—	—
Min. TR(ms)	—	2423	1691	1691	—	—	—
Min. TR/TI(ms)	—	—	—	—	—	—	—
BW in EPI freq. dir. (Hz)	—	—	—	—	—	—	—
SPIR offset act./default (Hz)	—	—	—	—	—	—	—
Local torso SAR	< 15 %	< 82 %	< 75 %	< 75 %	< 12 %	< 12 %	< 12 %
Whole body SAR / level	< 0.5 W/kg / normal	< 2.6 W/kg / 1st level	< 2.4 W/kg / 1st level	< 2.4 W/kg / 1st level	< 0.4 W/kg / normal	< 0.4 W/kg / normal	< 0.4 W/kg / normal
SED	0.0 kJ/kg	< 0.5 kJ/kg	< 0.3 kJ/kg	< 0.4 kJ/kg	< 0.1 kJ/kg	< 0.1 kJ/kg	0.0 kJ/kg
Coil Power	—	79 %	71 %	71 %	12 %	12 %	11 %
Max B1+rms	0.88 uT	2.07 uT	1.98 uT	1.98 uT	0.80 uT	0.80 uT	0.78 uT
PNS / level	78 % / normal	57 % / normal	55 % / normal	59 % / normal	75 % / normal	75 % / normal	65 % / normal
dB/dt	49.7 T/s	72.9 T/s	76.7 T/s	75.0 T/s	66.5 T/s	66.5 T/s	60.0 T/s
Sound Pressure Level (dB)	9.95438004	14.925498	10.8632116	11.2681313	14.8457556	14.8457556	13.8253288
Boil-off (hPa/h)	—	—	—	—	—	—	—
Max Boil-off (hPa/h)	—	—	—	—	—	—	—
Boil-off (hPa)	—	—	—	—	—	—	—

Disclaimer: The MR Protocols information is provided by independent physicians and technologists and is not intended to provide definitive or complete medical instructional or treatment information. Products referenced herein should be used consistently with their package inserts, and any devices referenced herein should be used consistently with the operation instructions provided by their respective manufactures. Northwest Imaging Forums, Inc. and Bracco Diagnostics Inc. do not offer medical advice and shall not be responsible for any person's reliance on this protocol.

Supported by an educational grant from Bracco Diagnostics Inc.

MR Gadopiclenol (VUEWAY®) Protocol for Philips MRI Scanners
Pediatric

IMAGING FORUMS, INC.

GEOMETRY

Scan Description	Survey	AX T2 FS MVXD	CO T2 FS MVXD	SG T2 FS MVXD	AX T1 FS mDIXON VANE 3D	AX T1 FS mDIXON VANE 3D POST CONTRAST	CO T1 FS 3D POST CONTRAST
Nucleus	—	—	—	—	—	—	—
SmartSelect	no	yes	yes	yes	yes	yes	yes
Coil 1 (exclude)	—	None	None	None	None	None	None
Uniformity	—	Synergy	CLEAR	CLEAR	CLEAR	CLEAR	CLEAR
FOV AP (mm)	—	200	135	200	180	180	120
FOV FH (mm)	300	199	200	200	156	156	180
FOV RL (mm)	300	200	200	189	180	180	180
stack RL (mm)	—	—	—	—	—	—	—
ACQ voxel size AP (mm)	—	1.10000002	—	1.29999995	1.29999995	1.29999995	1.29999995
ACQ voxel size FH (mm)	1.75781298	—	1.29999995	1.29999995	1.29999995	1.29999995	1.29999995
ACQ voxel size RL (mm)	3.515625	1.10000002	1.29999995	—	1.29999995	1.29999995	1.29999995
Slice thickness (mm)	15	4	3	4	—	—	—
Recon voxel size AP (mm)	—	0.833333313	—	0.892857134	1.02272725	1.02272725	1
Recon voxel size FH (mm)	1.5625	—	0.892857134	0.892857134	1	1	1.02272725
Recon voxel size RL (mm)	1.5625	0.833333313	0.892857134	—	1.02272725	1.02272725	1.02272725
Fold-over suppression	no	no	no	no	no	no	no
Slice oversampling	—	—	—	—	default	default	default
oversample factor	—	—	—	—	—	—	—
ENCASE enable	—	—	—	—	no	no	no
RF select. FOS	—	—	—	—	—	—	—
Reconstruction matrix	192	240	224	224	176	176	176
SENSE	no	yes	yes	yes	yes	yes	yes
Preduction (AP)	—	—	—	—	1	1	1
MB SENSE	—	—	—	—	—	—	—
CS-SENSE	no	—	—	—	—	—	—
reduction	—	—	—	—	—	—	—
denoising	—	—	—	—	—	—	—
k-t Acceleration	no	no	no	no	no	no	no
Stacks	3	1	1	1	1	1	1
current	—	—	—	—	—	—	—
type	parallel	parallel	parallel	parallel	—	—	—
slices	5	40	34	38	156	156	120
slice gap	user defined	user defined	user defined	user defined	—	—	—
gap (mm)	0	1	1	1	—	—	—
slice orientation	coronal	transverse	coronal	sagittal	transverse	transverse	coronal
fold-over direction	RL	—	—	—	—	—	—
fat shift direction	P	—	—	—	—	—	—
Minimum number of packages	—	1	1	1	—	—	—
Multi-chunk	—	—	—	—	no	no	no
Slice scan order	default	default	default	default	—	—	—
Stack scan order	ascend	—	—	—	—	—	—

Disclaimer: The MR Protocols information is provided by independent physicians and technologists and is not intended to provide definitive or complete medical instructional or treatment information. Products referenced herein should be used consistently with their package inserts, and any devices referenced herein should be used consistently with the operation instructions provided by their respective manufactures. Northwest Imaging Forums, Inc. and Bracco Diagnostics Inc. do not offer medical advice and shall not be responsible for any person's reliance on this protocol.

Supported by an educational grant from Bracco Diagnostics Inc.

Stack sort order	no	—	—	—	—	—	—
O-MAR	—	no	no	no	—	—	—
Large table movement	no	no	no	no	no	no	no
Plan Align	no	no	no	no	no	no	no
REST slabs	0	0	0	0	0	0	0
type	—	—	—	—	—	—	—
thickness (mm)	—	—	—	—	—	—	—
position	—	—	—	—	—	—	—
gap	—	—	—	—	—	—	—
power	—	—	2	2	—	—	—
Interactive positioning	no	no	no	no	no	no	no
Allow table movement	—	—	—	—	—	—	—
Patient position	head first	head first	head first	head first	head first	head first	head first
Patient body position	head first	head first	head first	head first	head first	head first	head first
Patient orientation	supine	supine	supine	supine	supine	supine	supine
Patient body orientation	supine	supine	supine	supine	supine	supine	supine

Disclaimer: The MR Protocols information is provided by independent physicians and technologists and is not intended to provide definitive or complete medical instructional or treatment information. Products referenced herein should be used consistently with their package inserts, and any devices referenced herein should be used consistently with the operation instructions provided by their respective manufactures. Northwest Imaging Forums, Inc. and Bracco Diagnostics Inc. do not offer medical advice and shall not be responsible for any person's reliance on this protocol.

Supported by an educational grant from Bracco Diagnostics Inc.

MR Gadopiclenol (VUEWAY®) Protocol for Philips MRI Scanners
Pediatric

IMAGING FORUMS, INC.

CONTRAST

Scan Description	Survey	AX T2 FS MVXD	CO T2 FS MVXD	SG T2 FS MVXD	AX T1 FS mDIXON VANE 3D	AX T1 FS mDIXON VANE 3D POST CONTRAST	CO T1 FS 3D POST CONTRAST
Scan type	Imaging	Imaging	Imaging	Imaging	Imaging	Imaging	Imaging
Scan mode	M2D	MS	MS	MS	3D	3D	3D
technique	FFE	SE	SE	SE	FFE	FFE	FFE
Modified SE	—	no	no	no	—	—	—
3D non-selective	—	—	—	—	—	—	—
loop order	—	—	—	—	—	—	—
Contrast enhancement	no	—	—	—	T1	T1	T1
Acquisition mode	cartesian	MultiVane	MultiVane	MultiVane	radial	radial	radial
Fast Imaging mode	TFE	TSE	TSE	TSE	TFE	TFE	TFE
3D VIEW	—	—	—	—	—	—	—
shot mode	single-shot	multishot	multishot	multishot	single-shot	single-shot	single-shot
TFE factor	85	—	—	—	96	96	74
TSE factor	—	50	28	28	—	—	—
startup echoes	default	0	0	0	default	default	default
shot interval	shortest	—	—	—	shortest	shortest	shortest
profile order	linear	linear	linear	linear	—	—	—
turbo direction	—	—	—	—	Z	Z	Z
DRIVE	—	yes	no	no	—	—	—
ultrashort	—	no	no	no	—	—	—
shift	—	—	—	—	—	—	—
fid reduction	—	default	default	default	—	—	—
Echoes	1	2	1	1	2	2	2
partial echo	no	no	no	no	no	no	no
shifted echo	no	—	—	—	no	no	no
TE	in-phase	shortest	user defined	user defined	shortest	shortest	shortest
(ms)	4.60384226	—	80	80	—	—	—
echospacing	—	—	—	—	—	—	—
(ms)	—	—	—	—	—	—	—
flyback	—	—	—	—	no	no	no
Flip angle (deg)	20	90	90	90	10	10	10
Refocusing control	—	constant	constant	constant	—	—	—
angle (deg)	—	80	80	80	—	—	—
bright fat reduction	—	no	no	no	—	—	—
TR	shortest	user defined	user defined	user defined	shortest	shortest	shortest
(ms)	—	2725	1860	1860	—	—	—
minimum (ms)	—	—	—	—	—	—	—
maximum (ms)	—	—	—	—	—	—	—
Halfscan	no	no	no	no	yes	yes	yes
factor	—	—	—	—	—	—	—
factor Y	—	—	—	—	—	—	—
factor Z	—	—	—	—	0.800000012	0.800000012	0.800000012

Disclaimer: The MR Protocols information is provided by independent physicians and technologists and is not intended to provide definitive or complete medical instructional or treatment information. Products referenced herein should be used consistently with their package inserts, and any devices referenced herein should be used consistently with the operation instructions provided by their respective manufactures. Northwest Imaging Forums, Inc. and Bracco Diagnostics Inc. do not offer medical advice and shall not be responsible for any person's reliance on this protocol.

Supported by an educational grant from Bracco Diagnostics Inc.

Water-fat shift	maximum	user defined	minimum	minimum	user defined	user defined	user defined
(pixels)	—	1	—	—	0.449999988	0.449999988	0.449999988
IR delay (ms)	—	—	—	—	—	—	—
dual	—	—	—	—	—	—	—
power	—	—	—	—	—	—	—
RF Shims	—	fixed	fixed	fixed	fixed	fixed	fixed
Shim	—	default	default	default	default	default	default
mDIXON	no	yes	no	no	yes	yes	yes
Fat suppression	no	no	SPAIR	SPAIR	no	no	no
strength	—	—	—	—	—	—	—
power	—	—	2	2	—	—	—
inversion delay	—	—	auto	auto	—	—	—
suppr. level	—	—	strong	strong	—	—	—
frequency offset	—	—	default	default	—	—	—
offset (Hz)	—	—	—	—	—	—	—
Grad Rev Fat suppr	—	no	no	no	—	—	—
Water suppression	no	no	no	no	no	no	no
TFE prepulse	no	—	—	—	no	no	no
slice selection	—	—	—	—	—	—	—
delay	—	—	—	—	—	—	—
(ms)	—	—	—	—	—	—	—
BB pulse	—	no	no	no	—	—	—
MTC	no	no	no	no	no	no	no
T2prep	no	—	—	—	no	no	no
echo time(ms)	—	—	—	—	—	—	—
refocusing pulses	—	—	—	—	—	—	—
APT	—	no	no	no	—	—	—
Custom prepulse	no	no	no	no	no	no	no
MDME	no	no	no	no	no	no	no
Diffusion mode	no	no	no	no	no	no	no
sequence	—	—	—	—	—	—	—
gradient expert mode	—	—	—	—	—	—	—
gradient overplus	—	—	—	—	—	—	—
direction	—	FH	—	—	FH	FH	FH
nr of b-factors	—	—	—	—	—	—	—
b-factor order	—	—	—	—	—	—	—
max b-factor	—	—	—	—	—	—	—
average high b	—	—	—	—	—	—	—
T1 mapping	no	no	no	no	no	no	no
Transmit channels	—	both	both	both	both	both	both
SAR mode	high	high	high	high	high	high	high
B1 mode	default	default	default	default	default	default	default
SAR allow first level	yes	yes	yes	yes	yes	yes	yes

Disclaimer: The MR Protocols information is provided by independent physicians and technologists and is not intended to provide definitive or complete medical instructional or treatment information. Products referenced herein should be used consistently with their package inserts, and any devices referenced herein should be used consistently with the operation instructions provided by their respective manufactures. Northwest Imaging Forums, Inc. and Bracco Diagnostics Inc. do not offer medical advice and shall not be responsible for any person's reliance on this protocol.

Supported by an educational grant from Bracco Diagnostics Inc.

Patient pregnancy	no	no	no	no	no	no	no
Patient WB SAR [W/kg]	0	0	0	0	0	0	0
Patient Head SAR [W/kg]	0	0	0	0	0	0	0
Patient max. dB/dt [T/s]	0	0	0	0	0	0	0
Max slewrate [T/m/s]	0	0	0	0	0	0	0
Max. B1+rms [uT]	0	0	0	0	0	0	0
PNS mode	moderate	high	high	high	high	high	high
Gradient mode	default	maximum	maximum	maximum	maximum	maximum	maximum
SofTone mode	no	no	no	no	no	no	no

MOTION

Scan Description	Survey	AX T2 FS MVXD	CO T2 FS MVXD	SG T2 FS MVXD	AX T1 FS mDIXON VANE 3D	AX T1 FS mDIXON VANE 3D POST CONTRAST	CO T1 FS 3D POST CONTRAST
Cardiac synchronization	no	no	no	no	no	no	no
Heart rate > 250 bpm	no	no	no	no	no	no	no
Respiratory compensation	no	no	trigger	trigger	no	no	no
Navigator respiratory comp	no	trigger and track	no	no	gate	gate	gate
type	—	parallel	parallel	parallel	lead	lead	lead
Flow compensation	yes	no	no	no	no	no	no
Temporal slice spacing	—	default	default	default	—	—	—
Motion smoothing	no	no	no	no	no	no	no
fMRI echo stabilisation	no	—	—	—	no	no	no
NSA	2	1	1	1	1	1	1
MRE enable	no	no	no	no	no	no	no

DYNAMIC ANGIO

Scan Description	Survey	AX T2 FS MVXD	CO T2 FS MVXD	SG T2 FS MVXD	AX T1 FS mDIXON VANE 3D	AX T1 FS mDIXON VANE 3D POST CONTRAST	CO T1 FS 3D POST CONTRAST
Angio / Contrast enh.	no	—	—	—	no	no	no
Quantitative flow	no	—	—	—	no	no	no
CENTRA	—	—	—	—	—	—	—
Manual start	no	no	no	no	no	no	no
Dynamic study	no	no	no	no	no	no	no
Arterial Spin labeling	no	no	no	no	no	no	no

Disclaimer: The MR Protocols information is provided by independent physicians and technologists and is not intended to provide definitive or complete medical instructional or treatment information. Products referenced herein should be used consistently with their package inserts, and any devices referenced herein should be used consistently with the operation instructions provided by their respective manufactures. Northwest Imaging Forums, Inc. and Bracco Diagnostics Inc. do not offer medical advice and shall not be responsible for any person's reliance on this protocol.

Supported by an educational grant from Bracco Diagnostics Inc.

MR Gadopiclenol (VUEWAY®) Protocol for Philips MRI Scanners
Pediatric

IMAGING FORUMS, INC.

POST PROCESSING

Scan Description	Survey	AX T2 FS MVXD	CO T2 FS MVXD	SG T2 FS MVXD	AX T1 FS mDIXON VANE 3D	AX T1 FS mDIXON VANE 3D POST CONTRAST	CO T1 FS 3D POST CONTRAST
Preparation phases	—	auto	auto	auto	full	full	full
Interactive F0	no	no	no	no	no	no	no
Quick Survey	default	default	default	default	default	default	default
SmartPlan survey	no	no	no	no	no	no	no
B0 field map	no	no	no	no	no	no	no
B1 field map	no	no	no	no	no	no	no
MIP/MPR	no	no	no	no	no	no	no
SWIp	no	—	—	—	no	no	no
Images	no	—	no	no	—	—	—
Autoview image	M	—	M	M	—	—	—
Calculated images	no	—	no	no	—	—	—
Reference tissue	Lung	Liver	Grey matter	Grey matter	Liver	Liver	Liver
Recon compression	—	No	No	No	Default	Default	Default
Preset window contrast	soft	soft	soft	soft	no	no	no
Reconstruction mode	immediate	immediate	immediate	immediate	immediate	immediate	immediate
Save raw data	no	no	no	no	no	no	no
Hardcopy protocol	no	no	no	no	no	no	no
Image filter	system default	system default	system default	system default	system default	system default	system default
Uniformity correction	no	no	no	no	yes	yes	yes
Geometry correction	default	default	default	default	none	none	none
EPI geometry correction	—	—	—	—	—	—	—
ADC correction	—	—	—	—	—	—	—
Elliptical k-space shutter	—	—	—	—	—	—	—
Free rotatable	no	no	no	no	no	no	no

Disclaimer: The MR Protocols information is provided by independent physicians and technologists and is not intended to provide definitive or complete medical instructional or treatment information. Products referenced herein should be used consistently with their package inserts, and any devices referenced herein should be used consistently with the operation instructions provided by their respective manufactures. Northwest Imaging Forums, Inc. and Bracco Diagnostics Inc. do not offer medical advice and shall not be responsible for any person's reliance on this protocol.

Supported by an educational grant from Bracco Diagnostics Inc.