NW
IMAGING FORUMS, INC.

PEDIATRIC BRAIN/ORBIT NON-RB 3T

INSTRUCTIONS AND PATIENT PREPARATION

Instructions	Non-RB orbit imaging (all other orbital lesions and tumors including optic nerve glioma) with whole-brain coverage.
Position	Head first supine
Coil	Head 20 or Head 64
Patient Preparation	Landmark per department pediatric protocol — pending sign-off.
Coverage	**Whole Brain coverage** — Pre-contrast Sag 3D T1 MPRAGE, TRAN rDWI; Post-contrast TRAN T2 FLAIR, SAG 3D T1 SPACE FS **Orbit coverage** — TRAN T2 TSE DIXON, COR T2 TSE DIXON, PRE- and POST TRAN T1 VIBE DIXON, and POST COR T1 TSE DIXON • Superior — Superior margin of midbrain • Inferior — Mid maxilla • Anterior — Tip of nose • Posterior — Mid pons • Left to Right — Skin to skin

CONTRAST AGENT PARAMETERS

Contrast Agent	Dose	Injection Rate	Dose Timing Method and Scan Delay
Gadopiclenol (VUEWAY)	0.05 mmol/kg	Hand Bolus or Power Injection	Per department pediatric contrast protocol — pending technologist/radiologist sign-off.
Comments	• Use minimal dose for patients with a GFR <60 mL/sec • Report agent administered and dose		

IMAGING PARAMETERS

Scan Description	SAG 3D T1 MPRAGE	TRANS rDWI	Tran _T2_TSE_Dixon DR	Tran_T1_VIBE_Dixon	Cor _T2_TSE_Dixon DR	Post Tran_T1_VIBE_Dixon	Post COR T1 TSE DIXON DR	Post TRAN T2 FLAIR DR	Post SAG 3D T1 SPACE FS
Sequence Name	tfl	resolve	tseR	fl	tseR	fl	tse	tir	spcR
Mode	Single Shot	Interleaved	Interleaved	—	Interleaved	—	Interleaved	Interleaved	—
PAT Mode	GRAPPA	—	GRAPPA	None	GRAPPA	None	GRAPPA	GRAPPA	CS
Dimension	3D	2D	—	3D	—	3D	—	2D	3D
Acceleration Factor PE	2	2	2	—	2	—	2	2	—
Acceleration Factor 3D	1	—	—	—	—	—	—	—	—

SCAN PARAMETERS

Scan Description	SAG 3D T1 MPRAGE	TRANS rDWI	Tran _T2_TSE_Dixon DR	Tran_T1_VIBE_Dixon	Cor _T2_TSE_Dixon DR	Post Tran_T1_VIBE_Dixon	Post COR T1 TSE DIXON DR	Post TRAN T2 FLAIR DR	Post SAG 3D T1 SPACE FS
TR	1540.0 ms	5160.0 ms	3220.0 ms	5.69 ms	4020.0 ms	5.69 ms	634.0 ms	10000.0 ms	900.0 ms
TE	2.76 ms	58 ms / 100 ms	83.00 ms	2.46 ms / 3.69 ms	83.00 ms	2.46 ms / 3.69 ms	12.00 ms	106.00 ms	19.00 ms
TI	900 ms	—	—	—	—	—	—	2600 ms	—
Flip Angle	15 deg	180 deg	160 deg	11 deg	160 deg	11 deg	160 deg	130 deg	—
Turbo Factor	132	—	21	—	21	—	5	12	44
Bandwidth	150 Hz/Px	875 Hz/Px	381 Hz/Px	740 Hz/Px	381 Hz/Px	740 Hz/Px	319 Hz/Px	349 Hz/Px	434 Hz/Px
EPI Factor	—	101	—	—	—	—	—	—	—
Fat Saturation Mode	—	Strong	Strong	—	Strong	—	Strong	Strong	Strong
Fat Suppression	Standard	Fat Saturation	Dixon	Dixon	Dixon	Dixon	Dixon	Fat Saturation	Fat Saturation
Gradient Mode	Normal	Fast	Normal	Fast	Normal	Fast	Normal	Fast	Fast
RF Pulse Type	Fast	Normal	Normal	Fast	Normal	Fast	Normal	Low SAR	Normal
Asymmetric Echo	Allowed	—	Off	Weak	Off	Weak	Off	—	—

Disclaimer: The MR Protocols information is provided by independent physicians and technologists and is not intended to provide definitive or complete medical instructional or treatment information. Products referenced herein should be used consistently with their package inserts, and any devices referenced herein should be used consistently with the operation instructions provided by their respective manufactures. Northwest Imaging Forums, Inc. and Bracco Diagnostics Inc. do not offer medical advice and shall not be responsible for any person's reliance on this protocol.

Supported by an educational grant from Bracco Diagnostics Inc.

MR Gadopiclenol (VUEWAY®) Protocol for Siemens MRI Scanners
Pediatric

IMAGING FORUMS, INC.

ACQUISITION PARAMETERS

Scan Description	SAG 3D T1 MPRAGE	TRANS rDWI	Tran _T2_TSE_Dixon DR	Tran_T1_VIBE_Dixon	Cor _T2_TSE_Dixon DR	Post Tran_T1_VIBE_Dixon	Post COR T1 TSE DIXON DR	Post TRAN T2 FLAIR DR	Post SAG 3D T1 SPACE FS
Base Resolution	192	168	320	160	320	160	320	256	256
Phase Resolution	100 %	100 %	75 %	90 %	75 %	90 %	70 %	80 %	100 %
Slice Resolution	75 %	—	—	75 %	—	75 %	—	—	75 %
Averages	1	—	1	1	1	1	1	1	1.0
Concatenations	1	1	3	1	3	1	5	2	1
b-value	—	0 s/mm²/800 s/mm²	—	—	—	—	—	—	—
FoV Read	192 mm	220 mm	150 mm	130 mm	150 mm	130 mm	150 mm	210 mm	192 mm
FoV Phase	100.0 %	100.0 %	106.3 %	117.5 %	100.0 %	117.5 %	100.0 %	93.8 %	100.0 %
Phase Encoding Dir.	A >> P	A >> P	R >> L	R >> L	R >> L	R >> L	R >> L	R >> L	A >> P
Phase Oversampling	20 %	20 %	10 %	15 %	10 %	15 %	20 %	0%	10 %
Slice Oversampling	10.0 %	—	—	100.0 %	—	100.0 %	—	—	16.7 %
Phase Partial Fourier	Off	Off	Off	Off	Off	Off	Off	Off	Off
Slice Partial Fourier	Off	—	—	Off	—	Off	—	—	Off
Flow Compensation	None	—	None	—	None	—	None	None	None

SCANNING RANGE

Scan Description	SAG 3D T1 MPRAGE	TRANS rDWI	Tran _T2_TSE_Dixon DR	Tran_T1_VIBE_Dixon	Cor _T2_TSE_Dixon DR	Post Tran_T1_VIBE_Dixon	Post COR T1 TSE DIXON DR	Post TRAN T2 FLAIR DR	Post SAG 3D T1 SPACE FS
Elliptical Filter	On	—	Off	Off	Off	Off	Off	Off	Off
Normalize	Prescan	Prescan	Prescan	Prescan	Prescan	Prescan	Prescan	Prescan	Prescan
Distortion Correction	2D	2D	2D	2D	2D	2D	2D	2D	2D
Slice Thickness	1.00 mm	4.0 mm	3.0 mm	0.8 mm	3.0 mm	0.8 mm	3.0 mm	4.0 mm	0.80 mm
Distance Factor	0%	0%	0%	0%	0%	0%	0%	0%	—
Slices	160	37	35	120	45	120	45	41	192
Acquisition Time	3:17 min	3:23 min	2:46 min	3:20 min	3:02 min	3:20 min	3:15 min	3:22 min	4:32 min

Disclaimer: The MR Protocols information is provided by independent physicians and technologists and is not intended to provide definitive or complete medical instructional or treatment information. Products referenced herein should be used consistently with their package inserts, and any devices referenced herein should be used consistently with the operation instructions provided by their respective manufactures. Northwest Imaging Forums, Inc. and Bracco Diagnostics Inc. do not offer medical advice and shall not be responsible for any person's reliance on this protocol.

Supported by an educational grant from Bracco Diagnostics Inc.